Product data sheet 5290-3 CT EAN-No. 4000896228904



Length: 291 mm · Torque range min-max: 10 – 60 Nm · Accuracy: 3 % · 📼



- Plug connection for holding insert tools
- Accuracy ± 3 % of the scale value (in the direction of actuation) Can be used for right-hand and left-hand threads by turning the insert tool
- Safe:
 - Haptic (close-gap release)
- Audible ("click" sound)
- Designed for rough workshop use
- Wide range of applications for controlled screw tightening in tyre service, bicycle, passenger car and commercial vehicle applications as well as various applications in trade and industry
- · Optimised sealing ring to protect against debris
- Handle with better grip for easier power transmission
- Adjustment aid with detent points for optimised operator guidance ensures safe and fast adjustment of the desired torque value by turning the handle
- Safe locking of the setting values by latching function on the rotary assembly
- Rotary assembly in triangular form prevents uncontrolled rolling
- Lock symbols indicate the respective locking status
- Can be attached to the cable loop through openings on the locking rotary assembly
- Easy-to-read, high-contrast scale
- Permanent readability due to laser marking on the scale bush
- With calibration certificate and serial number
- Delivery in sturdy Hexa plastic box
- DIN EN ISO 6789-2:2017
- Made in Germany
- Output: 💷
- Dimensions / length: 291 mm
- Net weight: 0.38 kg
- For manual operation
- Torque range min-max: 10-60 Nm
- Accuracy: 3 %
- Scale Graduation (Nm): 0.5 Nm
- Effective torque wrench length min.-max. I1 mm: 275.10 261.00 Nm
- Actuating force min.-max. kg: 3.7-23.4 kg
- Inside micrometre of the adjusting tool I2: 30 mm
- Head Thickness a: 17 mm
- Head Width b: 22 mm
- No more turning back: Continuous long-term tests on the 5000 series and 6000 series of HAZET torque wrenches have demonstrated that turning them back to the smallest scale value is not required





SYSTEM 5200-3 CT insert torque wrench \cdot 2.5 – 200 Nm \cdot with insert square \cdot accuracy \pm 3 %

Order Number

5290-3 CT